

SOP: SP022

Preparation of Alditol Acetate Derivatives

Materials and Reagents:

1. Sample, 20 to 50 µg
2. Rhamose standard, 10 mg/ml in B & J water
3. Fucose standard, 10 mg/ml in B & J water
4. Ribose standard, 10 mg/ml in B & J water
5. Arabinose standard, 10 mg/ml in B & J water
6. Xylose standard, 10 mg/ml in B & J water
7. Mannose standard, 10 mg/ml in B & J water
8. Galactose standard, 10 mg/ml in B & J water
9. Glucose standard, 10 mg/ml in B & J water
10. Myo-inositol standard, 10 mg/ml in B & J water
11. Scyllo-inositol standard, 1 mg/ml in B & J water
12. Trifluoroacetic acid, concentrated, 1 ml ampule
13. Water, Burdick & Jackson HPLC-grade
14. Methanol, Burdick & Jackson HPLC-grade
15. Sodium borodeuteride (NaBD₄), solid
16. Ammonium hydroxide, concentrated
17. Ethanol, absolute
18. Acetic acid, glacial
19. Acetic anhydride, 2 ml ampule
20. Chloroform, Burdick & Jackson HPLC-grade
21. Savant speed-vac
22. 13 x 100 mm glass tubes (twice as many as number of samples plus standard)
23. 13 mm PTFE-lined lids (as many as number of samples plus standard)
24. Capillary pipettor, 0-100 µl
25. Glass capillary pipets, 100 µl
26. Heat block, 120°C
27. Air bath
28. Glass Pasteur pipets
29. Rubber Pasteur pipet bulb
30. Dessicator
31. Vortex
32. Benchtop centrifuge
33. Capillary pipettor, 0-10 µl
34. Glass capillary pipets, 10 µl

Protocol:

1. ____ Transfer each sample into a cleaned 13 x 100 mm glass tube.
2. ____ Combine 25 µg of each sugar standard, except scyllo-inositol, in a new 13 x 100 mm glass tube.
3. ____ Completely dry all samples and the neutral sugar standard on the savant (note 1).
4. ____ Add 250 µl of 2M TFA to each sample and the standard (notes 2 and 3).
5. ____ Cap each tube tightly and place in 120°C heat block for two hours.
6. ____ Remove samples from the heat block and let cool to room temperature.
7. ____ Add 10 µg of scyllo-inositol to each sample and the standard (note 4).
8. ____ Completely dry the contents of each tube on the air bath (note 5).

9. ____ Add approximately 100 μ l of methanol to each sample (note 6).
10. ____ Completely dry on air bath.
11. ____ Repeat steps 9 and 10.
12. ____ Make NaBD₄ solution and add 250 μ l to each sample (note 7).
13. ____ Cap each tube and let sit on benchtop overnight (note 8).
14. ____ Add two drops of glacial acetic acid to each sample (note 9).
15. ____ Add 200 μ l of 10% acetic acid in methanol solution to each sample (note 10).
16. ____ Dry completely on air bath.
17. ____ Repeat steps 15 and 16.
18. ____ Add approximately 100 μ l of methanol to each sample (note 6).
19. ____ Completely dry on air bath.
20. ____ Repeat steps 18 and 19.
21. ____ Add 100 μ l of acetic anhydride from ampules to each sample.
22. ____ Cap each tube and heat at 120°C for two hours in a heating block.
23. ____ Remove samples from the heat block and let cool to room temperature.
24. ____ Completely dry on air bath.
25. ____ Add 1 ml of water to each sample.
26. ____ Add 2 ml of chloroform to each sample.
27. ____ Cap each tube and mix by vortexing vigorously.
28. ____ Centrifuge at 2,500 x g, 4°C for five minutes.
29. ____ Transfer the lower, organic layer from each sample into new 13 x 100 mm glass tubes and discard water layer (note 11).
30. ____ Completely dry on air bath.
31. ____ Sample is now ready for GC analysis.

Notes:

1. See SOP SP005 for use of savant
2. Use a fresh ampule of concentrated trifluoroacetic acid (TFA) to make 2M TFA for each alditol acetate preparation. Concentrated TFA is 12.98 M; use 153 μ l of acid to every 847 μ l of water to make each 1 ml of 2M TFA. Make in a glass container.
3. Use only glass capillary pipets to transfer liquids from this step on, as the GC will detect plastic components from Pipetman tips.

4. Scyllo-inositol is used as an internal standard for the GC to calculate the amount of neutral sugar in the sample.
5. See Air Bath SOP SP031.
6. This approximation is five drops from a glass Pasteur pipet with a rubber bulb. It is not necessary to add a specified amount of methanol, only enough to saturate the sample and remove any residual water from the sample.
7. The NaBD₄ solution is 10 mg of NaBD₄ in 1 ml of 1 M NH₄OH in ethanol, and must be freshly made prior to use. NaBD₄ is extremely hygroscopic and must be kept in a dessicator. 1 M NH₄OH in ethanol must be made fresh every two months. To make, add 6.6 ml of concentrated NH₄OH to 93.4 ml of absolute ethanol; as NH₄OH is caustic, make in a chemical fume hood.
8. The reduction reaction is complete in one hour, but overnight reduction provides the best results.
9. The addition of glacial acetic acid should cause the sample to bubble and fizz, indicating the required excess of reducing agent was present.
10. Make 10 % acetic acid in methanol as any other standard v/v solution; glacial acetic acid is caustic, so make in chemical fume hood.
11. It is important not to contaminate the organic layer with debris from the water layer. To do this, expel several drops of air while passing the Pasteur pipet through the water layer until the tip is in the organic layer. It is better to leave a small amount of organic layer in the tube than to risk water contamination.

Reference:

McNeil, M., D. Chatterjee, S. W. Hunter, and P. J. Brennan. 1989. Mycobacterial glycolipids: isolation, structures antigenicity and synthesis of neoantigens. *Methods Enzymology*. 179: 215-242.