

## SOP: SP044

### Sequence Alignment by Vector NTI Contig Express

#### Materials and Reagents:

1. Dell Dimension XPS T450 Rm C210
2. Electronic sequence chromatograms from forward, reverse, and any internal primers (note 1)
3. Electronic sequence of gene of interest in recombinant plasmid map (see SOP SP043) or from Tuberculist at <http://genolist.pasteur.fr/TubercuList/>, or other sources (eg GenBank).
4. Folder containing all sequence files and serving as save site for Contig(s) to be generated by this protocol

#### Protocol:

1. \_\_\_\_ Create a folder (note 2) containing electronic text and chromatograms compiled for sequence of interest.
2. \_\_\_\_ For pMRLB plasmid sequences, route folder to Dell desktop → Contract Related folder → rec plasmid info folder → specific plasmid folder.
3. \_\_\_\_ For sequences not (yet) assigned TB Contract reagent status, route folder to personal domain on T drive or other server.
4. \_\_\_\_ Open Vector NTI 9 – design window opens.
5. \_\_\_\_ Pull down Assemble (in toolbar) → select Contig Express → Open New Assembly Project. Contig Express Project window opens.
6. \_\_\_\_ Pull down Project (in toolbar) → select Add Fragments → select ABI files for sequence of interest.
7. \_\_\_\_ Migrate to folder created for this sequencing project in step 1.
8. \_\_\_\_ Open "All Files". ABI chromatogram list will appear.
9. \_\_\_\_ Highlight and Open ABI chromatogram file(s) in contig project window.
10. \_\_\_\_ In Vector NTI design window, pull down File → Local Database. Exploring window opens.
11. \_\_\_\_ Select subset with recombinant plasmid map file.
12. \_\_\_\_ Highlight plasmid map file and drag to contig project window.
13. \_\_\_\_ Alternative (steps 14 - 23): Create Vector NTI sequence file from Tuberculist (or other sequence source.)
14. \_\_\_\_ Download Tuberculist DNA sequence ± upstream/downstream sequence.
15. \_\_\_\_ Open As Windows default.
16. \_\_\_\_ Save as seq file, eg. Rv1908c.seq to personal domain.
17. \_\_\_\_ In Vector NTI Exploring window, pull down DNA/RNA in top toolbar.
18. \_\_\_\_ Press Import → Molecule from Text File.
19. \_\_\_\_ Press FastA – OK.

20. \_\_\_\_ Navigate to file. Select file. Open.
21. \_\_\_\_ Insert in DNA/RNA Main or in desired subset. New File window appears.
22. \_\_\_\_ Fill in gene name in General tab window. eg Rv1908c.seq
23. \_\_\_\_ Click OK. This file is available to join the contig.
24. \_\_\_\_ Drag file from project subset window to Contig Project window.
25. \_\_\_\_ In Contig Express Project window, highlight sequences to be aligned using Shift Click.
26. \_\_\_\_ In toolbar, press Assemble Selected (icon with arrow over three bars). Contig 1 is created.
27. \_\_\_\_ Blue arrows by sequences indicates contig has been assembled.
28. \_\_\_\_ Grey arrow indicates no match.
29. \_\_\_\_ Select Contig 1 icon to display alignment.
30. \_\_\_\_ Left window shows alignment fragments. Click Contig 1 to add description, change name.
31. \_\_\_\_ Right window shows graphic of file overlaps as bars. Green line shows region of poor alignment.
32. \_\_\_\_ Lower window shows base by base alignment in text.
33. \_\_\_\_ Clicking file name highlights that file in both graphic and text view.
34. \_\_\_\_ Contig 1, the consensus sequence, is shown in lower panel.
35. \_\_\_\_ Dots indicate unreadable sequence.
36. \_\_\_\_ Red sequence indicates mismatch.
37. \_\_\_\_ Add viewable sequence chromatograms to resolve mismatches (steps 38 – 41).
38. \_\_\_\_ Click in text pane. Pull down View → Show All Chromatograms. Chromatogram(s) are displayed below the corresponding sequence.
39. \_\_\_\_ View rating for one nucleotide by holding cursor over peak. Trace values appear.
40. \_\_\_\_ If value is high for a base reported as an “n”, change it in the text line.
41. \_\_\_\_ Note true mismatches (lack of agreement in sequence from two primers).
42. \_\_\_\_ Submit plasmid for resequencing of region(s) of unresolved mismatch.
43. \_\_\_\_ Information in all panes (description, graphic, text) can be printed if desired.

**Notes:**

1. Sequence is acquired from vendors (eg. SeqWright, Macromolecular Resources) as text and electronic chromatograms. In this protocol, electronic chromatograms are required for alignment with recombinant plasmid map (created in SOP SP043 Recombinant Plasmid Map Design – Vector NTI) or with gene sequence from Tuberculist database or other sources.

2. Label the folder "Contig.construct name.date" eg: Contig, Rv.1908,pET23b. 2-23-05. Store the contig folder in personal domain on T drive (or other server) unless it is added to the data assembly for TB contract plasmids. See SOP RM006 pMRLB Plasmid Documentation and SOP RM007 Non-pMRLB Plasmid Documentation for detailed listing of documentation requirements.

**References:**

Vector NTI 9 Online Help