

SOP: SP061.2

Modified: 6/16/2014 by GD

In-Solution Digestion of Proteins

Materials and Reagents:

1. ___ Protein of interest
2. ___ 0.65 ml Deplasticized Eppendorf tubes (note 1)
3. ___ Denaturing solution (6 M GuanHCl in 0.2M Tris-HCl, pH 8.6)
4. ___ Vortexer
5. ___ Reducing solution (0.4 M DTT, prepared in denaturing solution) (Note 2)
6. ___ Alkylating solution (18 mg iodoacetamide in 500 μ L 10mM ambic) (Note 2)
7. ___ 37°C incubator
8. ___ 10mM ammonium bicarbonate
9. ___ Slide-A-Lyzer MINI Dialysis Unit (Thermo #69552)
10. ___ Dialysis tank
11. ___ Savant Speed-Vac
12. ___ Modified trypsin sequencing grade (Roche Molecular Biochemicals, catalog # 1418 025) (Note 2)
13. ___ Acetonitrile (MS grade)
14. ___ Burdick and Jackson Water (note 3)
15. ___ Buffer A: 96.9% water, 3% ACN and 0.1% formic acid
16. ___ Autosampler vials with lids (Note 10)

Protocol:

Day 1:

1. ___ Savant Speed-Vac 100 μ g of the protein of interest to minimal volume (2-5 μ l) in a deplasticized Eppendorf tube. (Note 3 & 4)
2. ___ Resuspend protein in 100 μ L of denaturing solution.
3. ___ Add 1 μ L of reducing solution, vortex briefly and incubate at room temperature for 3 hours.
4. ___ Add 10 μ L of alkylating solution and incubate in the dark at 37°C for 30 minutes.
5. ___ Transfer samples to Slide-A-Lyzer unit and dialyze samples using 10 mM ammonium bicarbonate. Incubate overnight at 4 °C.

Day 2:

6. ___ Remove samples from dialysis units and rinse membrane with additional 10 mM ambic, adding total sample to a new deplasticized Eppendorf tube.
7. ___ Dry sample in Savant Speed-Vac down to minimal volume (2-5 μ l).
8. ___ Dissolve 25 μ g of trypsin in 300 μ l of 0.2 M ammonium bicarbonate (notes 5, 6 & 7).
9. ___ Add trypsin solution to the sample at a ratio of 50:1 (sample:trypsin), for 100 μ g of protein add 2 μ g (24 μ L) of trypsin
10. ___ Add acetonitrile to a final concentration of 10%
11. ___ Incubate the digest reaction overnight at 37°C.

Day 3:

12. ____ Terminate the reaction by adding 1/10th the digest volume of 10% TFA.
13. ____ Savant Speed-Vac and dry (note 8) (note 11).
14. ____ Samples **must** be cleaned prior to submission for LC-MS analysis. Appropriate clean-up methods may include C-18 zip-tips or spin columns.
15. ____ After cleaning savant Speed-Vac and dry the samples.
16. ____ Add buffer A to the sample and mix well. Sample should be at final concentration of 500 ng/μL - 1 μg/μL for the LTQ (note 9).
17. ____ Centrifuge the sample for 5 minutes.
18. ____ Pipet the sample out of the tube, being sure not to touch the bottom of the tube with the pipet tip (there will be a small volume left in the tube).
19. ____ Pipet sample into an autosampler vial, being careful not to introduce any bubbles into the vial.
20. ____ Place the cap on the autosampler vial and label the vial (note 10).
21. ____ Store the samples at -20°C until analysis by LC-MS-MS.

Notes:

1. __ Eppendorf tubes are deplasticized by filling with 60% acetonitrile-0.1% TFA, followed by mixing, incubation at room temperature for 1 hour and decanting of the solution. This process is repeated two times for each tube and the tubes are dried in the Savant Speed-Vac.
2. __ These solutions should be made fresh and discarded after use.
3. __ This procedure can be used with less than 100 μg of protein, adjust reagents accordingly.
4. __ See SOP SP005 for use of the savant.t.
5. __ This procedure can be used with proteases other than trypsin, however the buffer for the digestions may differ for other proteases.
6. __ All solutions should be made using Burdick and Jackson water.
7. __ When preparing a new vial of trypsin, be sure to date and initial it. The reconstituted trypsin should be stored at 4°C and can be used for up to two weeks.
8. __ Do not allow to dry completely. Leave 1-2 μl of liquid in the bottom of the tube.
9. __ Ideal final concentration varies greatly for the TQ-S depending on sample purity and purpose.
10. __ Do not use the blue auto sampler vial lids; use the orange ones
11. __ To do a double trypsin digestion, repeat steps 8-12 after using Savant Speed-Vac to dry sample; to do a single trypsin digestion, proceed to step 14.

References:

Hellman, U., C. Wernstedt, J. Gopez, and C. H. Heldin. 1995. Improvement of an "In-Gel" digestion procedure for the micro preparation of internal protein fragments for amino acid sequencing. *Anal Biochem* 224:451-5.

Rosenfeld, J. *et al.* 1992. In gel digestion of proteins for internal sequence analysis after 1 or 2-D gel electrophoresis. *Anal Biochem* 202: 173-179.

