

## SOP: M013

### Mycobactin HR Media and Plates

#### Materials

MilliQ water	1 L disposable filter (Nalgene)
MgCl <sub>2</sub>	Magnetic stir bar and plate
10 M NaOH	2 x 2 L bottles
L-asparagine	UltraPure agarose (Invitrogen)
2 M H <sub>2</sub> SO <sub>4</sub>	Large culture plates (15 ml)
ZnSO <sub>4</sub>	Class II Biosafety hood
MnCl <sub>2</sub>	Ziplock bags
KH <sub>2</sub> PO <sub>4</sub>	100 ml bottle

#### Protocol

1. \_\_\_\_ Place 950 ml MilliQ water in 2 L bottle reserved for mycobactin work. Add 20 ml glycerol and stir until homogenous.
2. \_\_\_\_ Adjust pH to 7.5 with concentrated NaOH, then add 5 g L-asparagine, 5 g KH<sub>2</sub>PO<sub>4</sub>.
3. \_\_\_\_ Adjust pH to 6.8 with H<sub>2</sub>SO<sub>4</sub> and add Zn<sup>2+</sup>, Mn<sup>2+</sup>, and Mg<sup>2+</sup> to obtain final concentrations shown below (Note 1).  
  

ZnSO <sub>4</sub>	0.46 µg/ml
MnCl <sub>2</sub>	0.10 µg/ml
MgCl <sub>2</sub>	0.04 mg/ml
4. \_\_\_\_ Allow solution to stir at least 30 min, filter into disposable container, then transfer to fresh 2 L bottle.
5. \_\_\_\_ Transfer 80 ml media to 100 ml bottle for autoclaving. This solution will be transported to BSL-3 facility for growth of Mtb.
6. \_\_\_\_ Add 18.4 g UltraPure agarose to 920 ml media in 2 L bottle. Autoclave at least 30 min.
7. \_\_\_\_ Working in Class II BSC, transfer warm media to large plates, 125 ml each. Once dry, transfer plates to ziplock bags for storage at 4°C.

#### Notes

- (1) This should make the ions more soluble prior to autoclaving media. Since these are such minute quantities, it is more efficient to aliquot from stock solutions made in endotoxin-free H<sub>2</sub>O.

#### References

RM Hall and C Ratledge (1982) A simple method for the production of mycobactin, the lipid-soluble siderophore from mycobacteria FEMS Microbiology Letters 15: 133-136.

