

**SOP: PP024.8**

**Modified: 2/27/2017 MCL**

**Purification of PhoS1 (38kDa protein) from the CFP**

**Materials and Reagents:**

1. Culture filtrate proteins (CFP) from *M. tuberculosis* or 50% cut supernatant (see SOP PP035)
2. Ammonium sulfate
3. Endotoxin-free water
4. Dialysis buffer (10mM ammonium bicarbonate, 1mM DTT)
5. Concanavalin A-Sepharose 4B resin (Sigma, C9017)
6. ConA Binding Buffer A: 50mM KH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl, 1mM each of MgCl<sub>2</sub>, CaCl<sub>2</sub>, MnCl<sub>2</sub> and DTT (note 1)
7. ConA Elution Buffer B: Binding buffer with 0.5 M Methyl D-mannose (note 2)
8. ConA Elution Buffer C: Binding buffer with 0.5 M Methyl  $\alpha$ -D-mannopyranoside (note 2)
9. 5 ml HiTrap phenyl sepharose column
10. Phen Seph Buffer A: 1M ammonium sulfate, 50mM Sodium Phosphate, 1mM DTT, pH 7.4 (note 3)
11. Phen Seph Buffer B: 50mM Sodium Phosphate, 1mM DTT, pH 7.4 (note 4)
12. 12 x 75 mm disposable culture tubes
13. Polypropylene oak ridge tubes
14. Dialysis tubing, 3500 MWCO
15. Dialysis tank
16. Biorad Econopump chromatography system
17. Sorvall centrifuge

**Protocol:**

- 1.\_\_\_\_\_ If starting from CFP, thaw at 4°C overnight.
- 2.\_\_\_\_\_ Add a stir bar, place on a stir plate and begin stirring. Slowly add ammonium sulfate to 40% saturation (note 5).
- 3.\_\_\_\_\_ Stir at room temperature until ammonium sulfate is completely dissolved, then transfer to 4°C for and stir for 4-16 hr.
- 4.\_\_\_\_\_ Remove stir bar and centrifuge the CFP/ammonium sulfate solution at 27,000 x g, 4°C for 1 hour.
- 5.\_\_\_\_\_ Decant to supernatant into a clean bottle and save the pellet at 4°C (note 6).
- 6.\_\_\_\_\_ While stirring, slowly add ammonium sulfate to 70% saturation (note 5), being sure to calculate based on the starting concentration of 40%.
- 7.\_\_\_\_\_ Stir at room temperature until ammonium sulfate is completely dissolved, then transfer to 4°C for and stir for 4-16 hr.
- 8.\_\_\_\_\_ Remove stir bar and centrifuge the CFP/ammonium sulfate solution at 27,000 x g, 4°C for 1 hour.
- 9.\_\_\_\_\_ Decant supernatant. The pellet is the 70% ammonium sulfate cut.
- 10.\_\_\_\_\_ Resuspend the pellet in 25-30 mL dialysis buffer and transfer to dialysis tubing.
- 11.\_\_\_\_\_ Dialyze at 4°C for 4-16 hours.
- 12.\_\_\_\_\_ Change dialysis buffer (7 L) and dialyze at 4°C for 4-16 hours.
- 13.\_\_\_\_\_ Change dialysis buffer (7 L) of 10 mM ammonium bicarbonate (without DTT) and dialyze at 4°C for 4-16 hours.

- 14.\_\_\_\_\_ Collect the protein solution from the dialysis tubing and rinse the tubing with a minimal volume of fresh 10 mM ammonium bicarbonate. Place the protein solution along with the rinse in a clean sample tube.
- 15.\_\_\_\_\_ Quantitate protein by BCA (see SOP: SP003), and lyophilize.
- 16.\_\_\_\_\_ Into an open column, pour a volume of Con-A resin which will give a ratio of 1 ml packed resin: 2 mg protein (note 7).
- 17.\_\_\_\_\_ Pack and equilibrate the column with 3-5 column volumes (CV) ConA Binding Buffer at a flow rate of 2 ml/min using the Biorad Econopump (note 8).
- 18.\_\_\_\_\_ Resuspend the dried sample in ConA Binding Buffer at a concentration of 1 mg/ml.
- 19.\_\_\_\_\_ Load the sample onto the column at 1 ml/min and pass the flow through back over the column. Collect as "ConA Flow Thru" (note 9).
- 20.\_\_\_\_\_ Run the following gradient at 1.5 ml/min:
  - 3 CV Binding Buffer A Wash (note 9)
  - 4 CV 0→45% Elution Buffer B (0.225M D-mannose)
  - 1 CV 45% Elution Buffer B hold
  - 2 CV 100% Elution Buffer B hold (0.5M D-mannose)
  - 2 CV 100% Elution Buffer C clean-up (0.5M Methyl  $\alpha$ -D-mannopyranoside)Collect the Wash as "ConA Wash", then collect 80 fractions starting at the gradient through the clean-up.
- 21.\_\_\_\_\_ Run 10  $\mu$ l of each fraction on a gel.
- 22.\_\_\_\_\_ Once it is determined which fractions contain 38kDa, run 10 ml of each in a western blot probed against Rv1411c (a common contaminant that is difficult to remove), being sure to include a positive control (CFP or recombinant Rv1411c).
- 23.\_\_\_\_\_ Pool all of the fractions that are predominantly 38kDa with little-to-no Rv1411c contamination.
- 24.\_\_\_\_\_ Concentrate the pool using an amicon ultra-15 and wash three times with 10 mM ammonium bicarbonate.
- 25.\_\_\_\_\_ Quantitate by BCA, visualize the pool by SDS-PAGE, and lyophilize.
- 26.\_\_\_\_\_ Suspend lyophilized ConA pool in Phenyl Sepharose Buffer A at a concentration of 1 mg/ml and stir at 4°C overnight (note 10).
- 27.\_\_\_\_\_ Pump 25 ml (5 CV) of endotoxin-free water through the HiTrap column at 1 ml/min to elute the storage buffer.
- 28.\_\_\_\_\_ Pump 25 ml of Buffer A through the column at 1 ml/min to equilibrate the column.
- 29.\_\_\_\_\_ Centrifuge the resuspended sample at 3000 rpm, 4°C, for 15 minutes to remove any precipitate.
- 30.\_\_\_\_\_ Pump the sample (supernatant) onto the column at a flow rate of 1 ml/min, then pass the flow through back over the column. Collect as "Phen Seph Flow Thru".
- 31.\_\_\_\_\_ Run the following gradient:
  - 2 CV Buffer A Wash (10 min)
  - 20 CV A→B Gradient (100 min)
  - 10 CV Buffer B Clean-Up (50 min)Collect the wash as "Phen Seph Wash", then collect 40 x 2.5 min fractions during the gradient. Collect the final clean up as "Phen Seph Clean-Up" (note 11).

- 32.\_\_\_\_\_ Run 10 µl of each fraction on a gel and check for Rv1411c contamination by western blot.
- 33.\_\_\_\_\_ Pool all fractions containing clean 38kDa.
- 34.\_\_\_\_\_ Concentrate using amicon ultra-15 centrifugal device and wash three times with 10 mM ammonium bicarbonate.
- 35.\_\_\_\_\_ Run BCA, gel, and blot using IT-23 (+) and Rv1411c (-) antibodies, for QC.
- 36.\_\_\_\_\_ Make aliquots (default quantity = 0.25 mg), lyophilize and store at -80°C.

**Notes:**

1. Con A Binding Buffer A:

To 80 ml of endotoxin-free water stirring on a stirplate, add the following:

KH <sub>2</sub> PO <sub>4</sub>	0.69 g
NaCl	2.92 g
MgCl <sub>2</sub> ·6H <sub>2</sub> O	20.3 mg
CaCl <sub>2</sub> ·2H <sub>2</sub> O	14.7 mg
MnCl <sub>2</sub> ·4H <sub>2</sub> O	19.8 mg
DTT	15.4 mg

After all reagents have gone into solution, add NaOH dropwise until the pH is 5.7. Transfer to a graduated cylinder and bring to final volume of 100 ml with endotoxin-free water.

2. Con A Elution Buffer B: Dissolve all the above reagents in 70 ml of endotoxin-free water plus D-mannose 9.01 g. Titrate to pH = 5.7 and bring to final volume of 100 ml.
- Con A Elution Buffer C: Dissolve all the above reagents in 70 ml of endotoxin-free water plus Methyl α-D-mannopyranoside 9.71 g. Titrate to pH = 5.7 and bring to final volume of 100 ml.

3. Buffer A:

To 70 ml of endotoxin-free water stirring on a stir plate, add the following reagents:

ammonium sulfate	13.22 g
NaH <sub>2</sub> PO <sub>4</sub>	0.114 g
Na <sub>2</sub> HPO <sub>4</sub>	1.0865 g
DTT	15.3 mg

When all reagents have gone into solution, titrate to pH = 7.4. Bring final volume to 100 ml with endotoxin-free water.

4. Buffer B:

To 70 ml of endotoxin-free water which is stirring on a stir plate, add the following reagents:

NaH <sub>2</sub> PO <sub>4</sub>	0.114 g
Na <sub>2</sub> HPO <sub>4</sub>	1.0865 g
DTT	15.3 mg

When all reagents have gone into solution, titrate to pH=7.4. Bring final volume to 100 ml with endotoxin-free water.

5. Determine the appropriate amount of ammonium sulfate using the calculator at <http://www.encorbio.com/protocols/AM-SO4.htm>.
6. The pellet is the 40% ammonium sulfate cut and is used as starting material for Ag85 (see SOP PP020 and PP021) and Mpt32 (SOP PP022). This material should be resuspended and dialyzed along with the 70% cut, and then freeze-dried and stored at -80°C for future use.
7. Measure approximately 1.5 ml resin slurry for every 1 ml of desired packed resin. For large preparations of CFP, it may be necessary to run more than one column, due to the amount of time necessary to run large columns.
8. The column can be packed and equilibrated the day before use, as long as there is sodium azide present in the buffer. All buffers containing sodium azide must be disposed of as hazardous waste.
9. Label as "ConA Flow Thru/Binding Buffer Wash – GroES Starting Material". Dialyze into 10 mM ammonium bicarbonate, lyophilize, and store at -80°C.
10. Even if the protein appears to go into solution it will often crash out of solution on the phenyl sepharose column, therefore overnight incubation is essential.

11. To clean the column for storage, pump 5 column volumes of water, then 5 column volumes of 20% ethanol. Cap both ends of the column and store at 4°C.